

PRODUCT INFORMATION

Common Name	SRF 388;SRF-388;SRF388
Synonyms	p28;IL30;IL-27;IL27A;IL-27A;IL27p28
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL27
Uniprot ID	Q8NEV9
Description	PE-conjugated Anti-IL27(casdozokitug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

