

PRODUCT INFORMATION

Common Name	BIIB-019;BIIB019;Dacliximab;RO 247375;RO247375;RO-247375;Zenapax;Zinbryta
Synonyms	p55;CD25;IL2R;IMD41;TCGFR;IDDM10;IL2-RA;IL-2-RA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL2RA
Uniprot ID	P01589
Description	PE-conjugated Anti-IL2RA(daclizumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

