

**PRODUCT INFORMATION**

<b>Common Name</b>	RG 6292;RG-6292;RG6962;RO 7296682;RO-7296682;RO7296682
<b>Synonyms</b>	p55;CD25;IL2R;IMD41;TCGFR;IDDM10;IL2-RA;IL-2-RA
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	IL2RA
<b>Uniprot ID</b>	P01589
<b>Description</b>	PE-conjugated Anti-IL2RA(vopikitug biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

DIMABIO CONFIDENTIAL

