

PRODUCT INFORMATION

Common Name	BI 655130;BI655130;BI-655130;SPEVIGO
Synonyms	FIL1;FIL1D;IL1F5;IL1L1;PSORP;IL1HY1;IL1RP3;IL36RA;IL-36Ra;PSORS14;FIL1(Delta)
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL36RN
Uniprot ID	Q9UBH0
Description	PE-conjugated Anti-IL36RN(spesolimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

