

PRODUCT INFORMATION

Common Name	BC005;BC-005;BC 005;MG-010;MG010;MG 010;MG-K10;MG K10
Synonyms	CD124;IL4RA;IL-4RA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	IL4RA
Uniprot ID	P24394
Description	PE-conjugated Anti-IL4RA(comekibart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

