

PRODUCT INFORMATION

Common Name	TBC
Synonyms	CD124;IL4RA;IL-4RA
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL4RA
Uniprot ID	P24394
Description	PE-conjugated Anti-IL4RA(telikibart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

