

PRODUCT INFORMATION

Common Name	Kevzara;REGN-88;REGN88;REGN 88;SAR-153191;SAR153191;SAR153191
Synonyms	IL6Q;gp80;CD126;HIES5;IL-6R;IL6RA;IL6RQ;IL-1Ra;IL-6RA;IL6QTL;IL-6R-1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL6R
Uniprot ID	P08887
Description	PE-conjugated Anti-IL6R(sarilumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

