

PRODUCT INFORMATION

Common Name	Effi-7, Effi7-h, MD-707
Synonyms	CD127; CDW127; IL-7R-alpha; IL7RA; ILRA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid □ PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG4
Reactivity	Human
Target	IL7RA
Uniprot ID	P16871
Description	PE-conjugated Anti-IL7RA(lusvertikimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

