

PRODUCT INFORMATION

Common Name	R2477;R-2477;REGN-2477;REGN2477
Synonyms	Inhibin beta A chain;INHBA;Activin A;EDF;FRP
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	INHBA
Uniprot ID	P08476
Description	PE-conjugated Anti-INHBA(garetosmab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

