

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | RZ-358;X 358;XOMA 247;XOMA-247;XOMA247;XOMA 358;XOMA-358;XOMA358;XPA.15.247 |
| Synonyms | HHF5;CD220 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Whole mAb |
| IgG type | IgG2 |
| Reactivity | Human |
| Target | INSR |
| Uniprot ID | P06213 |
| Description | PE-conjugated Anti-INSR(ersodetug biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

