

PRODUCT INFORMATION

Common Name	TBC
Synonyms	kidney associated antigen 1; kidney-associated antigen 1; MGC78738; RU2 antisense gene protein; RU2AS
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid □ PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	KAAG1
Uniprot ID	Q9UBP8
Description	PE-conjugated Anti-KAAG1(olintatug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

