

PRODUCT INFORMATION

Common Name	ABC-008,ABC008
Synonyms	2F1;CLEC15A;MAFA;MAFA-2F1;MAFA-L;MAFA-LIKE
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1(E239D, M241L)
Reactivity	Human
Target	KLRG1
Uniprot ID	Q96E93
Description	PE-conjugated Anti-KLRG1(ulviprubart biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

