

PRODUCT INFORMATION

Common Name	GSK-4074386;GSK4074386;TSR033;TSR033;TSR-033
Synonyms	LAG-3;CD223
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	LAG3
Uniprot ID	P18627
Description	PE-conjugated Anti-LAG3(encelimumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

