

PRODUCT INFORMATION

Common Name	BI-754111;BI754111;BI 754111
Synonyms	LAG-3;CD223
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	LAG3
Uniprot ID	P18627
Description	PE-conjugated Anti-LAG3(miptenalimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

