

PRODUCT INFORMATION

Common Name	AZD0171;AZD-0171;MSC-1;MSC1
Synonyms	CDF;DIA;HILDA;MLPLI
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	LIF
Uniprot ID	P15018
Description	PE-conjugated Anti-LIF(falbikitug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

