

PRODUCT INFORMATION

Common Name	AEVI-002;AEVI002;AEVI 002;AEVI-2;AEVI2;AEVI 2;AVTX-002;AVTX002;AVTX 002;AVTX-2;AVTX2;AVTX 2;CERC 002;CERC-002;CERC002;CERC 2;CERC-2;CERC2;MDGN-002;MDGN002;MDGN 002;MDGN 2;MDGN-2;MDGN2;SAR 252067;SAR252067;SAR-252067
Synonyms	LTg;CD258;HVEML;LIGHT;TNFSF14
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	LIGHT
Uniprot ID	O43557
Description	PE-conjugated Anti-LIGHT(quisovalimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

