

PRODUCT INFORMATION

Common Name	HZN-7734;HZN7734;MEDI-7734;MEDI7734;VIB-7734;VIB7734
Synonyms	ILT7;CD85g
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	LILRA4
Uniprot ID	P59901
Description	PE-conjugated Anti-LILRA4(daxdilimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

