

PRODUCT INFORMATION

Common Name	TBC
Synonyms	B4;ILT3;LIR5;CD85K;ILT-3;LIR-5;HM18
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	LILRB4
Uniprot ID	Q8NHJ6
Description	PE-conjugated Anti-LILRB4(rolistobart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

