

PRODUCT INFORMATION

Common Name	LA480, LY-2875358, LY2875358
Synonyms	MET
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG4(F114A,L115A)
Reactivity	Human
Target	MET
Uniprot ID	P08581
Description	PE-conjugated Anti-MET(emibetuzumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

