

PRODUCT INFORMATION

Clone ID	1F9
Target	Mouse IgG
Synonyms	N/A
Host Species	Goat
Description	PE-conjugated Anti-Mouse IgG antibody(1F9), Goat mAb
Delivery	Under Development
Uniprot ID	N/A
IgG type	Goat IgG
Clonality	Monoclonal
Reactivity	Mouse
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	N/A
Usage	Research use only
Conjugate	PE-conjugated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

DIMA BIO CONFIDENTIAL

