

## PRODUCT INFORMATION

|   |   |
|---|---|
| <b>Clone ID</b>                         | 7B12  |
| <b>Target</b>                           | Mouse IgG   |
| <b>Synonyms</b>                         | N/A   |
| <b>Host Species</b>                     | Goat  |
| <b>Description</b>                      | PE-conjugated Anti-Mouse IgG antibody(7B12),<br>Goat mAb          |
| <b>Delivery</b>                         | 3-4 weeks   |
| <b>Uniprot ID</b>                       | N/A   |
| <b>IgG type</b>                         | Goat IgG  |
| <b>Clonality</b>                        | Monoclonal  |
| <b>Reactivity</b>                       | Mouse   |
| <b>Applications</b>                     | Flow Cyt  |
| <b>Recommended Dilutions</b>            | Flow Cyt 1:100  |
| <b>Purification</b>                     | Purified from cell culture supernatant by affinity chromatography |
| <b>Formulation &amp; Reconstitution</b> | Liquid□PBS with 0.05% Proclin300, 1% BSA                          |
| <b>Storage &amp; Shipping</b>           | Store at 2°C-8°C for 6 months                                     |
| <b>Background</b>                       | N/A   |
| <b>Usage</b>                            | Research use only   |
| <b>Conjugate</b>                        | PE-conjugated   |

