

**PRODUCT INFORMATION**

<b>Common Name</b>	PankoMab,PankoMab-GEX,hPankoMab,DS-3939□DS-3939a
<b>Synonyms</b>	ADMCKD, ADMCKD1, CA 15-3, CD227, EMA, H23AG, KL-6, MAM6, MCD, MCKD, MCKD1, MUC-1, MUC-1/SEC, MUC-1/X, MUC1/ZD, PEM, PEMT, PUM
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	N/A
<b>Host Species</b>	Humanized
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	MUC1
<b>Uniprot ID</b>	P15941
<b>Description</b>	PE-conjugated Anti-MUC1(gatipotuzumab biosimilar) mAb
<b>Delivery</b>	3-4 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

