

PRODUCT INFORMATION

Common Name	PF-04383119;PF04383119;PF-4383119;PF4383119;RI-624;RI624;RI 624;RN-624;RN624;RN 624
Synonyms	NGFB;HSAN5;Beta-NGF;NGF;NGFB;-NGF
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG2
Reactivity	Human
Target	NGF
Uniprot ID	P01138
Description	PE-conjugated Anti-NGF(tanezumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

