

PRODUCT INFORMATION

Common Name	XMT-153;XMT153;XMT 153
Synonyms	PULAM;NPTIIb;NaPi2b;NAPI-3B;NAPI-IIb
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	NPT2B
Uniprot ID	O95436
Description	PE-conjugated Anti-NPT2B(upinitatug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

