

PRODUCT INFORMATION

Common Name	JS-002;JS002;JS 002;JS-2;JS2;JS 2
Synonyms	FH3;PC9;FHCL3;NARC1;LDLCQ1;NARC-1;HCHOLA3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	PCSK9
Uniprot ID	Q8NBP7
Description	PE-conjugated Anti-PCSK9(ongericimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

