

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | SHR-1209;SHR1209;SHR 1209 |
| Synonyms | FH3;PC9;FHCL3;NARC1;LDLCQ1;NARC-1;HCHOLA3 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | PCSK9 |
| Uniprot ID | Q8NBP7 |
| Description | PE-conjugated Anti-PCSK9(recaticimab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

