

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | AGEN-2034;AGEN2034;RebmAb-700;RebmAb700 |
| Synonyms | PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | N/A |
| Host Species | Whole mAb |
| IgG type | IgG4 |
| Reactivity | Human |
| Target | PD-1 |
| Uniprot ID | Q15116 |
| Description | PE-conjugated Anti-PD-1(balstilimab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

