

PRODUCT INFORMATION

Common Name	APL 501;APL-501;APL501;CBT 501;CBT501;CBT-501;GB 226;GB-226;GB226
Synonyms	PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	PD-1
Uniprot ID	Q15116
Description	PE-conjugated Anti-PD-1(geptanolimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

