

PRODUCT INFORMATION

Common Name	CT-011;CT011;CT 011;CT-11;CT11;CT 11;MDV 9300;MDV-9300;MDV9300
Synonyms	PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	PD-1
Uniprot ID	Q15116
Description	PE-conjugated Anti-PD-1(pidilizumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

