

PRODUCT INFORMATION

Common Name	CHS-007;CHS007;CHS 007;CHS-7;CHS7;CHS 7;JS001;JS001;JS-001;JS1;JS-1;JS 1;SO-001;SO001;SO01;SO1;SO-1;SO 1;TAB-001;TAB001;TAB001;TAB1;TAB-1;TAB 1;Tripleitriumab;Triprizumab;Tuoyi
Conjugate	PE-conjugated
Synonyms	PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	N/A
IgG type	IgG4
Reactivity	Human
Target	PD-1
Uniprot ID	Q15116
Description	PE-conjugated Anti-PD-1(toripalimab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

