

PRODUCT INFORMATION

Common Name	CK 301;CK301;CK-301;TG 1501;TG1501;TG-1501
Synonyms	B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	PD-L1
Uniprot ID	Q9NZQ7
Description	PE-conjugated Anti-PD-L1(cosibelimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

