

**PRODUCT INFORMATION**

<b>Common Name</b>	A-167;A 167;A167;HBM-9167;HBM9167;HBM 9167;KL-A167;KL A167
<b>Synonyms</b>	B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1
<b>Applications</b>	N/A
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	N/A
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	PD-L1
<b>Uniprot ID</b>	Q9NZQ7
<b>Description</b>	PE-conjugated Anti-PD-L1(tagitanlimab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

DIMABIO CONFIDENTIAL

