

PRODUCT INFORMATION

Common Name	PF-06523435,hu24
Synonyms	CCK4
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1(E239D, M241L)
Reactivity	Human
Target	PTK7
Uniprot ID	Q13308
Description	PE-conjugated Anti-PTK7(cofetuzumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

