

PRODUCT INFORMATION

Clone ID	37A11
Target	Rabbit IgG
Synonyms	N/A
Host Species	Rabbit
Description	PE-conjugated Anti-Rabbit IgG antibody(37A11), Goat mAb
Delivery	3-4 weeks
Uniprot ID	N/A
IgG type	Goat IgG
Clonality	Monoclonal
Reactivity	Rabbit
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	N/A
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

