

PRODUCT INFORMATION

Common Name	BA302
Synonyms	NTRKR2
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric Humanized
IgG type	IgG1(E239D, M241L)
Reactivity	Human
Target	ROR2
Uniprot ID	Q01974
Description	PE-conjugated Anti-ROR2(ozuriftamab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

