

PRODUCT INFORMATION

Common Name	SEG101, SelG1, crizanlizumab-tmca
Synonyms	P-Selectin;CD62P;SELP;GMP-140
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG2
Reactivity	Human
Target	SELP
Uniprot ID	P16109
Description	PE-conjugated Anti-SELP(crizanlizumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

