

PRODUCT INFORMATION

Common Name	SelK-2;SelK2
Synonyms	CLA;CD162;PSGL1;PSGL-1;Selectin P ligand
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG2
Reactivity	Human
Target	SELPLG
Uniprot ID	Q14242
Description	PE-conjugated Anti-SELPLG(cepeprubart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

