

## PRODUCT INFORMATION

<b>Common Name</b>	VX15/2503;VX 15;VX-15;VX 15;VX 15/2503;VX-15/2503
<b>Synonyms</b>	A8;GR3;BB18;CD100;COLL4;SEMAJ;coll-4;C9orf164;M-sema-G
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG4
<b>Reactivity</b>	Human
<b>Target</b>	SEMA4D
<b>Uniprot ID</b>	Q92854
<b>Description</b>	PE-conjugated Anti-SEMA4D(pepinemab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

DIMABIO CONFIDENTIAL

