

PRODUCT INFORMATION

Common Name	GS0189;GS-0189;GS189;GS-189;FSI189;FSI-189
Synonyms	BIT;MFR;P84;SIRP;MYD-1;SHPS1;CD172A;PTPNS1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	SIRP α
Uniprot ID	P78324
Description	PE-conjugated Anti-SIRP α (epacmarstobart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

