

PRODUCT INFORMATION

Common Name	LY-2541546;LY2541546;TST-002;TST002
Synonyms	CDD;VBCH;DAND6;SOST1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	SOST
Uniprot ID	Q9BQB4
Description	PE-conjugated Anti-SOST(blosozumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

