

**PRODUCT INFORMATION**

<b>Common Name</b>	BPS-804;BPS804;BPS 804;MOR-05813;MOR05813;MOR 05813;MOR-5813;MOR5813;MOR 5813
<b>Synonyms</b>	CDD;VBCH;DAND6;SOST1
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG2
<b>Reactivity</b>	Human
<b>Target</b>	SOST
<b>Uniprot ID</b>	Q9BQB4
<b>Description</b>	PE-conjugated Anti-SOST(setrusumab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

DIMABIO CONFIDENTIAL

