

PRODUCT INFORMATION

Common Name	XmAb-18087
Synonyms	SS-2-R;SS2-R;SS2R;SST2;SRIF-1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric Humanized
IgG type	IgG1
Reactivity	Human
Target	SSTR2
Uniprot ID	P30874
Description	PE-conjugated Anti-SSTR2(tidutamab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

