

PRODUCT INFORMATION

Common Name	MORAb-004
Synonyms	CD248; CD164L1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric Humanized
IgG type	IgG1(E239D,M241L,S286F)
Reactivity	Human
Target	TEM1
Uniprot ID	Q9HCU0
Description	PE-conjugated Anti-TEM1(ontuxizumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

