

PRODUCT INFORMATION

Common Name	SRK-181;SRK 181;SKR181
Synonyms	CED;LAP;DPD1;TGFB;IBDIMDE;TGFbeta;TGF-beta1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	TGFB1
Uniprot ID	P01137
Description	PE-conjugated Anti-TGFB1(linavonkibart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

