

PRODUCT INFORMATION

Common Name	BGB-A1217;BGB A1217;WCD-118;WCD118;WCD 118
Synonyms	VSIG9;VSTM3;WUCAM
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TIGIT
Uniprot ID	Q495A1
Description	PE-conjugated Anti-TIGIT(ociperlimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

