

PRODUCT INFORMATION

Common Name	BGB-A425;BGB A425
Synonyms	TIM3;CD366;KIM-3;SPTCL;TIMD3;Tim-3;TIMD-3;HAVcr-2
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TIM3
Uniprot ID	Q8TDQ0
Description	PE-conjugated Anti-TIM3(surzebiclimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

