

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | OPN-305 |
| Synonyms | CD282;TIL4 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | N/A |
| Host Species | Chimeric |
| IgG type | IgG4 |
| Reactivity | Human |
| Target | TLR2 |
| Uniprot ID | O60603 |
| Description | PE-conjugated Anti-TLR2(tomaralimab biosimilar) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

