

**PRODUCT INFORMATION**

<b>Common Name</b>	EB 05;EB05;EB-05;EB5;EB-5;EB5;EB-07;EB07;EB-07;EB-7;EB7;EB7;NI-0101;NI 0101;NI0101;NI-101;NI101;NI 101
<b>Synonyms</b>	TOLL;CD284;TLR-4;ARMD10
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	TLR4
<b>Uniprot ID</b>	O00206
<b>Description</b>	PE-conjugated Anti-TLR4(paridiprubart biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

