

PRODUCT INFORMATION

Common Name	AGX101
Synonyms	M3S1;TAAL6
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	TM4SF1
Uniprot ID	P30408
Description	PE-conjugated Anti-TM4SF1 (AGX101 biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

