

PRODUCT INFORMATION

Common Name	AGX101
Conjugate	PE-conjugated
Synonyms	M3S1;TAAL6
Applications	Flow Cyt
Recommended Dilutions	
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	TM4SF1
Uniprot ID	P30408
Description	PE-conjugated Anti-TM4SF1 (AGX101 biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.






Figure 1. Flow cytometry analysis with 1 μ l/test PE-conjugated Anti-TM4SF1 (AGX101 biosimilar) mAb on SNU-5 cell line (Blue histogram) or 293T cell line (Red histogram).

DIMABIO CONFIDENTIAL

